CONTROLLED DOCUMENT

WIRION™
Embofic Protection System
Instructions for Use

Excluding U.S.

Document No. P1-5-0721 Rev.11

22 March 2016

WIRION™ Embolic Protection System Instructions for Use

INTENDED USE

The WIRION Embolic Protection System is indicated for use as an embolic protection system to contain and remove embolic material during cardiovascular interventions.

CONTRAINDICATIONS

- Do not use laser devices with the System
- Any contraindication for PCI

PRODUCT DESCRIPTION

The WIRION Embolic Protection System is a temporary embolic protection system, filtering distal to the intervention site. It is provided sterile and for single use. The system is a rapid exchange, pre-loaded filter that can be used with commercially available 0.014" guide wires. After the guidewire of choice is properly located in the target vessel, the filter is delivered, locked and deployed on the wire, at any location according to physician's discretion. The locking mechanism is remotely activated, allowing the stand-alone filter unit to be positioned in any location along the vessel without compromising wire performance and support. At the end of the procedure, the filter is retrieved along with the wire.

Embolic material are blood clots, plaque or other solid substrates that break loose during the treatment procedure and enter into the blood stream resulting in downstream obstruction and damage to the end organ. Embolic protection devices collect and remove emboli and debris, reducing incidence of no reflow and damage to the distal tissue.

SYSTEM COMPONENTS

The system is supplied sterile (sterilized by ETO).

The system consists of a Delivery Catheter and a Retrieval Catheter;

Delivery Catheter: Contains and delivers a Filter Unit, locks it onto the guide wire using the Activating Handle and deploys it in the target vessel (Fig. 1). The Filter Unit is made of a filter membrane and a Nitinol frame. A 5cc syringe with a soft needle is supplied for flushing the system prior to use.

Fig. 1: WIRION™ Delivery Catheter

Retrieval Catheter: Retrieves the Filter Unit with the captured emboli and debris after the therapeutic procedure is completed.


Fig. 2: WIRION™ Retrieval Catheter

Quant.	Description
1	Delivery Catheter
1	Retrieval Catheter
1	5cc Flushing Syringe
1	Instructions for Use

Table 1: System Components

Device Characteristics	
Blood vessel diameter	3.5-6.0mm
Filter unit length	22.5 mm

Table 2: Device Characteristics

WARNINGS

Only physicians who have received appropriate training and are familiar with the principles, clinical applications, complications and hazards commonly associated with interventional procedures should use this system.

- Intended for single-use only. Do not reuse. Reuse of this device may lead to its failure and subsequent patient injury, illness or death.

- Do not resterilize as this can compromise device performance and increase the risk of cross contamination due to in appropriate processing.
- In cases the delivery system was retracted from the patient’s body, prior to locking the filter on the guide wire, it will not be re-used and will be discarded. The system will be replaced by a new one.
- The WIRION system is provided sterile. Do not use if the sterile barrier is damaged or if the package label is missing. If the barrier is compromised, contact your Gardia Medical representative.
- Refer to the instructions supplied with any interventional device to be used in conjunction with the System for their intended use, contraindications, warnings and precautions.
- The appropriate antiplatelet and anticoagulation therapy should be administered pre- and post-procedure to minimize the risk of embolism and thrombosis.
- Introduce and advance devices slowly to prevent air embolism or trauma to the vasculature.
- Do not use excessive force when attempting to cross the lesion with the System
- Do not use the product after the “Use By” date specified on the package

PRECAUTIONS

- Inspect the catheters Prior to use (e.g., bends, kinks, integrity). Do not use if damaged.
- Maintain the sterile contents of the pouch in the sterile field during the procedure.
- Confirm angiographically that the vessel diameter is appropriate for treatment with the System.
- Confirm angiographically that there is a sufficient vessel length in which to place the Filter Unit, and to maintain adequate distance between the Filter Unit proximal side and the stent delivery device or other interventional devices.
- Ensure use of a guiding catheter or introducer sheath to allow system support during filter deployment; allow no more than 35cm length outside of the guiding catheter/sheath.
- Use a 6F to 8F guiding catheter with a 0.07” (1.78mm) to 0.088” (2.24mm) inner diameter respectively, or a 6F introducer sheath with a 0.079” (2mm) minimum inner diameter,.
- Always use the System Retrieval Catheter to remove the Filter Unit with the guide wire. Pulling an open Filter Unit through a stent should be a last resort for retrieval as it may lead to entanglement of the filter with the stent.
- Always advance or retract the catheter slowly under fluoroscopic imaging. If resistance is felt and/or observed, determine the cause and take any necessary remedial action.
- Advancing or torquing the catheters against resistance could result in vessel trauma.
- Use caution when advancing or retracting the catheters through a deployed stent as this may cause stent entanglement or stent dislocation.
- If blood flow through the filter is slow, the filter may be removed and replaced. However, at the operator’s discretion, the procedure may be completed before removing the device.
- If there is no blood flow through the filter, remove the filter, using the Retrieval Catheter. If further intervention is required, replace with a new system.
- If intervention is required in additional vessels, use a new System.

POTENTIAL ADVERSE EVENTS

Possible adverse events include, but are not limited to the following:

- Angina
- Bleeding complications
- Bradycardia or arrhythmias, including ventricular fibrillation or tachycardia
- Congestive heart failure
- Damage to or dislocation of the implanted stent(s)
- Death
- Detachment and/or implantation of a component of the device
- Drug reaction, allergic reaction to contrast media, medications or device materials
- Emergent surgery
- Embolization of air, tissue, thrombus or other embolic debris
- End organ ischemia, vessel thrombosis or spasm
- Hypotension/hypertension
- Infection (local or systemic)
- Myocardial infarction
- No-reflow resulting from reduced blood flow through the Filter Unit
- Puncture site complications (i.e., vessel occlusion, hemorrhage, hematoma, pseudoaneurysm or arteriovenous fistula)
- Renal insufficiency, kidney failure, hematuria
- Stroke/cerebrovascular accident (CVA), transient ischemic attack (TIA) or seizure
- Vessel damage, dissection, occlusion, aneurysm, perforation, rupture or injury

PACKAGING & STORAGE

- Store the system in a dry place at room temperature.
- WIRION system is supplied in two pouches in a single box.

PHYSICIAN TRAINING

The System should be used by physicians thoroughly trained in percutaneous intravascular techniques and procedures. Training is required and may include in-service training and monitoring of initial case(s). These requirements will be based on individual physician experience.

INSTRUCTIONS FOR USE

1. Open the Delivery Catheter pouch using standard sterile handling procedures and place the catheter in the sterile field.
2. Take the Delivery Catheter out of the packaging coil.
3. Visually inspect the catheter for faults, bends or kinks. Special attention should be taken when inspecting the catheter distal part; the Filter Unit can now be seen inside the transparent amber tube, make sure that the soft tip is protruding out of the amber tube and the tube is intact.
4. Inspect the Activating Handle and make sure that the safety sticker is intact.

CAUTION: Do not use the Delivery Catheter if any of the above faults are found.

CAUTION: Do not peel away the safety sticker until the delivery and positioning of the catheter is completed

5. Remove the metal stylet placed in the catheter distal lumen.
6. Flush the Delivery Catheter with saline using the supplied 5cc flushing syringe.

CAUTION: Make sure the device is thoroughly flushed and liquid runs out of the Rx port and from the soft tip zone.

DELIVERY AND POSITIONING OF THE SYSTEM DELIVERY CATHETER

NOTE: At this stage of the procedure, a standard guide wire is already placed in the target vessel.

1. Carefully insert the Delivery Catheter onto the guide wire.
2. Advance the catheter on the guide wire until the guide wire exits the Rx port located about 38cm proximal to the catheter tip.
3. Under fluoroscopic imaging, advance the Delivery Catheter to the target vessel and place the Filter Unit distal to the lesion.
4. Once the Delivery Catheter crosses the lesion and is placed in the required position, slide a wire torque tool along the guide wire and secure it against the hemostasis valve.

NOTE: Position the proximal marker of the filter system more than 1cm distal to the expected position of the stent delivery system tip. Placement should be done according to the landing zone recommendations illustrated in Fig. 3.

NOTE: If extra landing area is available, it is suggested to place the Filter Unit distal to the minimum recommended distance.

NOTE: The Filter Unit and the marker bands are visible under fluoroscopic imaging.

NOTE: The technique of advancing the Delivery Catheter is similar to a standard PTCA balloon or stent device.

WARNING: Do not use excessive force when attempting to cross the lesion with the Delivery Catheter


Fig. 3: Landing zone clearance

LOCKING AND DEPLOYMENT OF THE DEVICE

NOTE: Once the Delivery Catheter is placed in the required position, the Filter Unit can be locked and deployed.


Fig. 4: Safety sticker removal


Fig. 5: Filter Marker Bands

1. Hold the Activating Handle and peel away the safety sticker (Fig 4).
2. Rotate the rotating section of the handle clockwise, until a red indicator becomes visible at the handle proximal end. The Filter Unit is now locked to the guide wire.
3. Deploy the Filter Unit by holding the guide wire in place with the guide wire torque tool pressed against the hemostasis valve while simultaneously retracting the Delivery Catheter. The three marker bands in the center of the Filter Unit will expand and indicate filter sealing position (Fig 5).
4. Continue retracting the Delivery Catheter, until the Rx port appears at the hemostasis port.
5. Remove the wire torque tool and carefully remove the distal segment of the Delivery Catheter from the wire while maintaining guide wire stability and Filter Unit position.

DURING THE INTERVENTION

1. Use the guide wire according to normal practice to track catheters and stent delivery devices to the target treatment site. Carefully maintain the position of the filter, particularly during device exchange.

WARNING: Always keep the open Filter Unit distal to the deployed stent. Pulling the Filter Unit into the stent area may lead to entanglement with the stent.

WARNING: Ensure that the guide wire is stabilized throughout the procedure. Failing to stabilize the guide wire could lead to inadvertent movement of the filter, resulting in guide wire entanglement and/or delay in the procedure.

2. Carefully remove all interventional devices, ensuring that the guide wire position is maintained during removal.

RETRIEVAL AND FILTER UNIT CAPTURE

CAUTION: When retrieving the Filter Unit, always use a WIRION Retrieval Catheter.

1. Open the Retrieval Catheter pouch using standard sterile handling procedures and place the catheter in the sterile field.
2. Take the Retrieval Catheter out of its packaging coil.
3. Remove the metal stylet placed in the catheter distal lumen by gently pulling its ring.
4. Flush the Retrieval Catheter with saline.
5. While maintaining the guide wire position distal to the lesion, advance the Retrieval Catheter over the guide wire past any deployed stent(s) until the tip of the Retrieval Catheter reaches the Filter Unit.
6. While holding the Retrieval catheter Handle (Fig. 2), clockwise rotate the Retrieving Knob by 90° and pull back (Fig. 6), until it reaches its final position (Fig. 7). The distal tip slides backwards as its marker band passes the collecting tube marker in the proximal direction. The catheter is now ready for retrieval.
7. Slide a wire torquer along the guide wire and secure it against the hemostasis valve.
8. Gently retrieve the Filter Unit into the Retrieval Catheter collecting tube. The Filter Unit and Retrieval Catheter tip will be visible under fluoroscopic imaging.
9. Carefully retract the Retrieval Catheter through the guide catheter or sheath up to the hemostasis valve.
10. Remove the Retrieval Catheter from the patient, ensuring that the hemostasis valve is sufficiently open prior to removal.

NOTE: According to physicians' decision, the Filter Unit can be fully or partially retrieved.

Fig 6: WIRION Retrieval Catheter- Handle starting position


Fig 7: WIRION Retrieval Catheter - retrieving position





CAUTION: Use caution when advancing or retracting the Retrieval Catheter through a deployed stent as this may cause Filter Unit/stent entanglement or stent dislocation.

WARNING: Do not pull excessively on the guide wire or the Retrieval Catheter.

NOTE: If any resistance is met during retraction of the guide wire and Filter Unit, slightly advance the guide wire and rotate the Retrieval Catheter before continuing to retract.

NOTE: If there is any resistance at the guiding catheter or sheath, retract the guide catheter or sheath and Retrieval Catheter together.

NOTE: After use, dispose of product and packaging in accordance with hospital administrative and/or local government policy

	Caution: See Instructions for Use
	Catalog Number
	Lot Number
	Sterilization using Ethylene Oxide
	Use By:
	EU Authorized Representative
	Do Not Reuse
	Warning: Do not use if package is damaged or opened
	Manufacturer
	Keep in cool dry location
	Caution: Federal (USA) law restricts this device to sale by or on order of a physician
	Blood Vessel Diameter
	Minimum and Maximum Outer Guiding Catheter Diameter
	System package sterile contents
	Guide wire insertion channel

Manufacturer:

Gardia Medical Ltd.
 2 Ha-Eshel St.
 Caesarea Industrial Park
 3088900, Israel
 Tel: + 972 - 4 - 6277166
 Fax: + 972 - 4 - 6277266

EC Authorized Representative:

CEpartner4U BV
 Esdoornlaan 13
 3951 DB Maarn
 The Netherlands
 Tel: +31- 343 -442524
 Fax: +31- 343 - 442162

Copyright 2012

WIRION® is a trademark of Gardia Medical Ltd.

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WIRION® jest znakiem towarowym firmy Gardia Medical Ltd.